Caramuel 38, 28011 Madrid • Tel. 91 464 94 50 - 91 464 36 00  
Fax. 91 464 62 58 • www.f-soria.es

FICHA TÉCNICA: 770495

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## REFERENCIAS BIBLIOGRAFICAS

1. Bauer, A.W., W.M.M. Kirby, J.C. Sherris, and M. Turck. 1966. Antibiotic susceptibility testing by a standardized single disk method. *Am. J. Clin. Pathol.* 45:493-496.
2. Ryan, K.J., F.D. Schoenknecht, and W.M.M. Kirby. 1970. Disc sensitivity testing. *Hospital Practice* 5:91-100.
3. Barry, A.L., F. Garcia, and L.D. Thrupp. 1970. An improved single-disk method for testing the antibiotic susceptibility of rapidly-growing pathogens. *Am. J. Clin. Pathol.* 53:149-158.
4. Mueller, J.H., and J. Hinton. 1941. A protein-free medium for primary isolation of the gonococcus and meningococcus. *Proc. Soc. Exp. Biol. Med.* 48:330-333.
5. Ericsson, H.M., and J.C. Sherris. 1971. Antibiotic sensitivity testing. Report of an international collaborative study. *Acta Pathol. Microbiol. Scand. Sec. B, Suppl.* 217.
6. National Committee for Clinical Laboratory Standards. 2000. Approved standard: M2-A7. Performance standards for antimicrobial disk susceptibility tests, 7<sup>th</sup> ed. National Committee for Clinical Laboratory Standards, Wayne, Pa.
7. Woods, G.L., and J.A. Washington. 1995. Antibacterial susceptibility tests: dilution and disk diffusion methods, p. 1327-1341. In P.R. Murray, E.J. Baron, M.A. Pfaller, F.C. Tenover, and R.C. Tenover (ed.), *Manual of clinical microbiology*, 6th ed. American Society for Microbiology, Washington, DC.
8. Thornsberry, C., T.L. Gavan, and E.H. Gerlach. 1977. Cumitech 6, New developments in antimicrobial agent susceptibility testing. Coordinating ed., J.C. Sherris. American Society of Microbiology, Washington, DC.
9. Koch, A.E., and J.J. Burchall. 1971. Reversal of the antimicrobial activity of trimethoprim by thymidine in commercially prepared media. *Appl. Microbiol.* 22:812-817.
10. Ferone, R., S.R.M. Bushby, J.J. Burchall, W.D. Moore, and D. Smith. 1975. Identification of Harper-Cawston factor as thymidine phosphorylase and removal from media of substances interfering with susceptibility testing to sulfonamides and diaminopyrimidines. *Antimicrob. Agents Chemother.* 7:91-98.
11. Reller, L.G., F.D. Schoenknecht, M.A. Kenny, and J.C. Sherris. 1974. Antibiotic susceptibility testing of *Pseudomonas aeruginosa*: selection of a control strain and criteria for magnesium and calcium content in media. *J. Infect. Dis.* 130:454-463.
12. Pollock, H.M., B.H. Minshew, M.A. Kenny, and F.D. Schoenknecht. 1978. Effect of different lots of Mueller-Hinton Agar on the interpretation of the gentamicin susceptibility of *Pseudomonas aeruginosa*. *Antimicrob. Agents Chemother.* 14:360-367.
13. D'Amato, R.F., and C. Thornsberry. 1979. Calcium and magnesium in Mueller-Hinton agar and their influence on disk diffusion susceptibility results. *Current Microbiol.* 2:135-138.
14. Wegner, D.L., C.R. Mathis, and T.R. Neblett. 1976. Direct method to determine the antibiotic susceptibility of rapidly growing blood pathogens. *Antimicrob. Agents Chemother.* 9:861-862.
15. Johnson, J.E., and J.A. Washington II. 1976. Comparison of direct and standardized antimicrobial susceptibility testing of positive blood cultures. *Antimicrob. Agents Chemother.* 10:211-214.

16. Waterworth, P.M., and M. Del Piano. 1976. Dependability of sensitivity tests in primary culture. *J. Clin. Pathol.* 29:179-184.
17. Hollick, G.E., and J.A. Washington II. 1976. Comparison of direct and standardized disk diffusion susceptibility testing of urine cultures. *Antimicrob. Agents Chemother.* 9:804-809.
18. Baker, C.N., C. Thornsberry, and R.W. Hawkinson. 1983. Inoculum standardization in antimicrobial susceptibility testing: evaluation of overnight agar cultures and the rapid inoculum standardization system. *J. Clin. Microbiol.* 17:450-457.
19. Hindler, J.A., and C.B. Anderbied. 1985. Effect of the source of Mueller-Hinton agar and resistance frequency on the detection of methicillin-resistant *Staphylococcus aureus*. *J. Clin. Microbiol.* 21:205-210.
20. Maskell, R., O.A. Okubadejo, R.H. Payne, and L. Pead. 1977. Human infections with thymine-requiring bacteria. *J. Med. Microbiol.* 11:33-45.
21. Haltiner, R.C., P.C. Migneault, and R.G. Robertson. 1980. Incidence of thymidine-dependent enterococci detected on Mueller-Hinton agar with low thymidine content. *Antimicrob. Agents Chemother.* 18:365-368.
22. Murray, B.E. 1990. The life and times of the Enterococcus. *Clin. Microbiol. Rev.* 3:46-65.
23. National Committee for Clinical Laboratory Standards. 1997. Approved standard: M7-A4 Methods for dilution antimicrobial susceptibility tests for bacteria that grow aerobically, 4th ed. National Committee for Clinical Laboratory Standards, Wayne, Pa.
24. Jorgensen, J.H., and J.D. Turnidge. 2003. Susceptibility test methods: dilution and disk diffusion methods. *In: Murray, P. R., E. J. Baron, J.H. Jorgensen, M. A. Pfaller, and R. H. Tenover (ed.). Manual of clinical microbiology, 8th ed. American Society for Microbiology, Washington, D.C.*

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