Caramuel 38, 28011 Madrid • Tel. 91 464 94 50 - 91 464 36 00  
Fax. 91 464 62 58 • www.f-soria.es

FICHA TÉCNICA: 770230

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Producto: **HAEMOPHILUS TEST MEDIUM AGAR**  
**PLACA DE 90 mm**

## USO

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## REFERENCIAS BIBLIOGRAFICAS

1. National Committee for Clinical Laboratory Standards. 2000. Approved standard: M2-A7. Performance standards for antimicrobial disk susceptibility tests, 7th ed. National Committee for Clinical Laboratory Standards, Wayne, Pa.
2. Thornsberry, C., T.L. Gavan, E.H. Gerlach. 1977. Cumitech 6, New developments in antimicrobial agent susceptibility testing. Coordinating ed., J.C. Sherris. American Society for Microbiology, Washington, D.C
3. Jorgensen, J.H., J.S. Redding, L.A. Maher, and A.W. Howell. 1987. Improved medium for antimicrobial susceptibility testing of *Haemophilus influenzae*. *J. Clin. Microbiol.* 25:2105-2113.
4. Jorgensen, J.H., A.W. Howell, and L.A. Maher. 1990. Antimicrobial susceptibility testing of less commonly isolated *Haemophilus* species using *Haemophilus* test medium. *J. Clin. Microbiol.* 28:985-988.
5. Baker, C.N., C. Thornsberry, and R.W. Hawkinson. 1983. Inoculum standardization in antimicrobial susceptibility testing: evaluation of overnight agar cultures and the rapid inoculum standardization system. *J. Clin. Microbiol.* 17:450-457.
6. Marsik, F., G. Evans, J. Fowler, and L. Thompson. 1989. Comparison of the BBL Prompt□ system, Abbott A-Just□, and visual method for the preparation of *Haemophilus influenzae* inoculum for the Bauer-Kirby procedure, abstr. C-67, p. 404. Abstr. 89th Annu. Meet. Am. Soc. Microbiol. 1989.
7. Doern, G.V., J.H. Jorgensen, C. Thornsberry, D.A. Preston, T.A. Tubert, J.S. Redding, and L.A. Maher. 1988. A national collaborative study of the prevalence of antimicrobial resistance among clinical isolates of *Haemophilus influenzae*. *Antimicrob. Agents Chemother.* 32:180-185.
8. Data on file. BD Diagnostic Systems, Sparks, MD, USA.
9. National Committee for Clinical Laboratory Standards. 1990. Approved standard: M2-A4. Performance standards for antimicrobial disk susceptibility tests, 4th ed. National Committee for Clinical Laboratory Standards, Villanova, Pa.

10. Azemun, P., T. Stull, M. Roberts, and A.L. Smith. 1981. Rapid detection of chloramphenicol resistance in *Haemophilus influenzae*. *Antimicrob. Agents Chemother.* 20:168-170.
11. Doern, G.V., G.S. Daum, and T.A. Tubert. 1987. In vitro chloramphenicol susceptibility testing of *Haemophilus influenzae*: disk diffusion procedures and assays for chloramphenicol acetyltransferase. *J. Clin. Microbiol.* 25:1453-1455.
12. National Committee for Clinical Laboratory Standards. 1984. Approved standard: M2-A3. Performance standards for antimicrobial disk susceptibility tests, 3rd ed. National Committee for Clinical Laboratory Standards, Villanova, Pa.
13. National Committee for Clinical Laboratory Standards. 1987. Second Informational Supplement: M100-S2. Performance standards for antimicrobial susceptibility testing. National Committee for Clinical Laboratory Standards, Villanova, Pa.
14. Doern, G.V. and R.N. Jones. 1988. Antimicrobial susceptibility testing of *Haemophilus influenzae*, *Branhamella catarrhalis*, and *Neisseria gonorrhoeae*. *Antimicrob. Agents Chemother.* 32:1747-1753.
15. Woods, G.L., and J.A. Washington. 1995. Antibacterial susceptibility tests: dilution and disk diffusion methods, p. 1327-1341. In P.R. Murray, E.J. Baron, M.A. Pfaller, F.C. Tenover, and R.H. Tenover (ed.), *Manual of clinical microbiology*, 6th ed. American Society for Microbiology, Washington, D.C. USA.

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